



CTM01

The **CTM01** Transport Container for Biological Samples is composed, in accordance with the OMS recommendations, of **two differentiated containers** that, concentrically, are housed one inside the other, plus a very resistant compact cardboard mandrel, that is housed in the inside the secondary container.

Specific for the transport of **Category A** samples.

This category, includes the set of substances with pathogenic agents, that can cause permanent damage or can endanger life, if exposure occurs.

Features

	OUTER PACKAGING	INNER PACKAGING
Material	Compact cardboard cylinder with an expanded polystyrene layer adapted to its inner wall.	Compact cardboard cylinder with embedded metal bottom.
Empty weight	73 g	35 g
Plastic cap weight	25 g	9 g
External diameter	99 mm	67 mm
Total height	170 mm (capless); 172 mm (with cap)	148 mm (capless); 149 mm (with cap)
Cobb1800	151,4 g/m ²	
Clasp type	Mobile lid on the upper part of the cylinder, made of rigid plastic and press-fitted to the body of the container, without a seal.	Mobile cover on the upper part of the cylinder, made of rigid plastic and snap-fitted to the cylinder body, without a seal.

Materials to pack

Container **CTM01** is indicated for the transport of all types of Infectious Biological Samples of those cataloged as Type 6.2.

Complies with Packing Instructions **P620** and **P650**. It can be used for both **CATEGORY A** (UN 2814) and **CATEGORY B** (UN 3373) infectious samples.

Sanitary Clinical Supplies

We are specialists in disposable sanitary supplies and equipment.

General description

Mandril	Made of compact cardboard with high resistance to compression.
Primary Containers	The packaging can contain an indeterminate number of primary containers depending on the type of container and size to be used (test tubes, flasks, etc).
Materials to pack	The CTM01 Container is indicated for the transport of all kinds of Infectious Biological Samples of those cataloged as Type 6.2. Complies with Packing Instructions P620 and P650 . It can be used for both CATEGORY A (UN 2814) and CATEGORY B (UN 3373) infectious samples.
Other considerations	BAR CODE: 8435104100103 CEE MARKING: Does not have. It is a logistic and non-sanitary product, therefore it does not require this marking.
Presentation	The CTM01 Containers are served in cardboard boxes of 50 units / box (minimum order).



Contact



(+34) 985 111 238



esdesla.com



info@esdesla.com



Asturias, Spain